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
Instruction Manual/Gebrauchsanleitung  
Mode d'emploi/Manual de instrucciones

*Leica*

# LEICA GALEN™ III Compound Microscope

## Instruction Manual

### Table of Contents

Section	Page
1.0 INTRODUCTION .....	2
2.0 UNPACKING .....	2
3.0 PREPARATION FOR USE .....	4
3.1 Objectives .....	4
3.2 Viewing Head .....	4
3.3 Condenser .....	4
3.4 Filter .....	4
3.5 Eyepieces .....	4
4.0 OPERATION .....	4
4.1 Auto Focus Stop .....	4
4.2 Setting Interpupillary Distance .....	6
4.3 Focusing .....	6
4.4 Blue Daylight Filter .....	8
4.5 Mechanical Stage .....	8
4.6 Operating Conditions .....	10
5.0 MAINTENANCE .....	10
5.1 Cleaning .....	10
5.2 Lamp Replacement .....	10
5.3 Service .....	10

### 1.0 INTRODUCTION

The LEICA GALEN™ III Compound Microscope is designed for use by clinical and education professionals as well as advanced science students, including those in veterinary and medical schools. A full range of accessories is available including adaptors for brightfield, darkfield, phase contrast, polarized light and photo- and videomicrography.

### 2.0 UNPACKING

Carefully remove the LEICA GALEN III from the packing carton. Before discarding any packaging, please check the following standard equipment list to assure that all components have been shipped.

1. Stand, including supporting arm, focusing mechanism, mechanical stage, nosepiece and mirror or integral illumination system
2. Viewing head (monocular, binocular or trinocular), as ordered
3. Objectives, as ordered
4. Eyepiece(s), as ordered
5. Abbe condenser with aperture diaphragm and swing-out filter holder (brightfield models)
6. Blue filter
7. Dust cover
8. Immersion oil

Optional accessories such as phase contrast kit, darkfield condenser, carrying case and camera accessories are shipped separately.

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## 3.0 PREPARATION FOR USE

### 3.1 Objectives

**CAUTION:** To prevent damage to the objectives during assembly to the nosepiece, rotate the coarse adjustment knob counter-clockwise for right hand adjustment to lower the stage.

### 3.2 Viewing Head

Remove the protective plastic cap from the viewing head. Loosen the viewing head clamp screw slightly and carefully seat the head onto the microscope arm. Lock the head into desired orientation by tightening the clamp screw.

### 3.3 Condenser

Loosen the knurled condenser clamp screw slightly and insert the condenser from below. Position the condenser so that the handle of the aperture diaphragm is conveniently accessible. Tighten the condenser clamp screw with your fingers. To eliminate light irregularity when using low power objectives such as 4X and 10X objectives, raise and lower the condenser using the condenser adjustment knob.

### 3.4 Filter

The filter holder is attached to the condenser allowing it to pivot into or out of the light path, as desired. Pivot the filter holder outward and insert the filter.

### 3.5 Eyepieces

Place the eyepiece(s) into the eyepiece tube(s). Each eyepiece has a built-in eyeguard that may be pulled out to reduce stray light or glare for eyeglass wearers.

The LEICA GALEN III is now assembled and ready for use.

## 4.0 OPERATION

Plug the microscope power cord into a suitable grounded electrical outlet. A grounded 3-wire cord is provided.

**CAUTION:** As a safety precaution, the power cord has been grounded to the microscope base. Never use an adapter between the power cord and the power source; it will render this feature ineffective.

Move the illuminator switch to the "on" position. To obtain the desired illumination, slide the light control rheostat to a suitable position. The degree of light intensity depends on individual conditions such as specimen contrast, objective magnification, individual eyesight, etc. Too little or too much light is not recommended. The energy saving illuminator indicator lamp remains lighted when the illuminator is in use.

### 4.1 Auto Focus Stop

The LEICA GALEN III is equipped with an auto focus stop to prevent the microscope stage from moving during certain applications. If the auto focus stop is required, after focusing tighten the auto focus knob by turning it clockwise. To assure free stage motion, turn the auto focus stop knob counterclockwise before focusing.

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## 4.2 Setting Interpupillary Distance

The LEICA GALEN III is designed to accommodate individual interpupillary distance (IPD). To assure that the IPD is set comfortably for individual use, grasp the microscope head with both hands and move the eyetubes together or apart until the full field of view is simultaneously visible through both eyepieces.

Rotate the diopter ring on the right eyepiece tube to the same setting. Place a specimen on the microscope stage. Close the left eye, or cover the left eyepiece with an opaque material and carefully focus the microscope with the fine adjustment knob. Close the right eye and without touching either adjustment knob, rotate the left diopter ring until the specimen image is in sharp focus for the left eye.

Adjusting the LEICA GALEN III for proper IPD is important for retaining objective parfocality as well as reducing operator fatigue.

### If the interocular and diopter ring adjustments are properly set:

- the specimen will have the same focus in the right and left eyepieces when either the coarse or fine adjustment knobs are used;
- the full field of view will be comfortably seen while your head is in a fixed position.
- only a fine adjustment will be necessary to bring an image into sharp focus as objectives are rotated for use.

Repeat the adjustments until all of the above results are realized. Once the instrument is properly set, note the reading on the interocular scale. This will assure the adjustments are set for your individual use each time.

**NOTE FOR EYEGLASS WEARERS:** *Optimum settings for the diopter rings may be slightly different for eyeglass wearers. If one eye requires greater correction, set the diopter ring on the eyepiece tube corresponding to your strong eye to the reading on the interocular scale. Close your weaker eye and carefully focus the microscope with the fine adjustment knob. Now close the strong eye and without touching either adjustment knob, rotate the diopter corresponding to the weak eye until the image is sharply focused.*

## 4.3 Focusing

Coarse Adjustment - Turn the right-hand side coarse adjustment knob counterclockwise. The stage will stop at the bottom of its travel. Turning the knob in the opposite direction will raise the stage.

Place a specimen slide on the center of the stage. Using a 10X objective and a 10X eyepiece for preliminary observation bring the specimen into focus with the coarse adjustment knob.

**CAUTION:** *To prevent damage to the specimen and objective:*

- Without using the eyepieces, directly look at the slide and raise the stage by slowly rotating the right-hand side coarse adjustment knob clockwise. Continue raising the stage until the cover glass nearly touches the objective or the stage reaches its upper limit.
- While viewing the specimen through the eyepieces, turn the coarse adjustment knob counterclockwise to lower the stage and bring the specimen into focus.

Adjust the Auto Focus Stop as described in Section 4.1. The other objectives will remain parfocal without danger of touching the specimen slide. If the stage moves too freely adjust the tension using the coarse adjustment tension knob by turning clockwise to tighten and counterclockwise to loosen it.

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#### 4.4 Blue Daylight Filter

If the light directed onto the specimen is excessively yellow, place the blue daylight filter in the condenser filter holder.

#### 4.5 Mechanical Stage

When a specimen is placed on the mechanical stage, the movement range is approximately 50mm x 75mm. The mechanical stage control knobs located directly to the right of the stage and in front of the coarse adjustment knob, should be rotated to change the viewing position of the specimen on the mechanical stage.

The scale located along the top, back section provides two readings: the main scale "m" measures in millimeters; the vernier scale is represented by "n".

$$1 \text{ division of } n = \frac{9 \text{ divisions of } m}{10}$$

**NOTE:** Figure 1 shows "n" positioned between 12mm and 13mm. The nearest point where "m" and "n" are aligned is used to determine the added millimeter division.

In this representation, Figure 1 shows 20 of m is aligned with 8 of n, or 0.8mm. Therefore, the reading is 12mm of m plus 0.8mm of n, or 12.8mm. This reading is the sum of the m indication, plus the reading at the closest aligned point.

**CAUTION:** Carefully and fully retract the mechanical stage when you have finished using the microscope. If the mechanical stage is fully extended either way, the stage could be damaged when attempting to put the microscope into the carrying case.

#### 4.6 Operating Conditions

Carefully following the preceding instructions, the LEICA GALEN III can be operated over a wide range of conditions with a minimum of adjustments. There are variables, however, which may be changed for use from specimen to specimen.

- The light level should be adjusted to the lowest level compatible with comfortable viewing. This will greatly extend the life of the lamp. The optimum lamp setting will change depending upon specimen density and objective magnification.
- When examining specimens of inherently low contrast, contrast can be improved by reducing the condenser numerical aperture (NA) to approximately 75% of the objective NA. This is best accomplished by removing one of the eyepieces and looking down the eyepiece tube at the illuminated back aperture of the objective. While observing this aperture, close down the condenser aperture diaphragm until the size of the objective back aperture is reduced to approximately 75% of its original size.
- Specimens with contrast so low that the specimen is invisible or barely visible often can be seen more clearly by closing down the condenser aperture diaphragm to a pinhole size.
- An alternate and often superior method of improving contrast is to use oblique illumination. This is accomplished by covering 1/2 or more of the condenser aperture with an opaque paper.


Figure 1 / Abb. 1 / Figure 1 / Figura 1

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Both of these methods of contrast enhancement introduce diffraction effects which reduce the ability of the microscope to resolve the fine details of the specimen. They are not recommended where a critical examination of the specimen is required. Instead, use phase contrast technique microscopy.

Accessories are available for converting a standard LEICA GALEN III for phase contrast microscopy. Contact your Leica representative for assistance with these accessories.

## 5.0 MAINTENANCE

Keeping the LEICA GALEN III clean, along with minimum routine maintenance will insure that the microscope has a long life.

### Important:

- Always protect the microscope with a dust cover when it is not in use.
- Always keep the objectives and eyepieces in place when storing the microscope.

## 5.1 Cleaning

The outer surface of the LEICA GALEN III should be cleaned from time to time with a soft damp cloth. Use only a mild soap solution to clean more stubborn dirt.

The outer surfaces of the optics should be inspected for dust and dirt. Use an air bulb to blow dust from the optics. If necessary, use a soft cloth or cotton swab dampened with a mild solution of soap and water to clean the recessed optical surfaces. Use alcohol or xylene only as a last resort.

**NOTE:** Avoid the use of solvents. They may cause problems with cemented optics. Solvents may also attract grease from the mounts making cleaning more difficult.

After each use of immersion oil, clean the objectives and stage with a soft cloth or lens tissue. Occasionally the underside of the mechanical stage may become coated with immersion oil and require cleaning.

Detach the mechanical stage by removing the three screws on the underside of the stage table. When cleaning oil from the underside of the stage be careful not to remove the stage slide paths or contaminate them with lubricant. Reinstall the mechanical stage with care to prevent damage to the slide path, rack and pinion.

## 5.2 Lamp Replacement

The only component in the LEICA GALEN III that may need periodic replacement is the illuminator lamp. The LEICA GALEN III utilizes a 6V, 20W tungsten halogen lamp (Part Number 1120).

**DANGER:** The glass envelope of the lamp may be extremely hot. Do not attempt to change the lamp before it is completely cooled or without wearing adequate skin protection.

Before replacing the lamp, unplug the instrument. Move the mechanical stage to the extreme right position and carefully place the microscope on its left side. Open the lamp door in the bottom plate by pulling the knob.

After the lamp has cooled, carefully remove it from its socket. Place the new lamp in its socket. Close the lamp door. Never operate the microscope illuminator unless the lamp door is securely closed.

## 5.3 Service

All optical and mechanical parts require periodic servicing. To keep the instrument performing properly and to compensate for normal wear, contact your local Leica representative to establish a regular schedule for servicing.

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Printed in USA ©1995, Leica Inc., Buffalo NY USA 10/95 317505-100 Rev. B
Leica Inc.